

## SPECTROTEST & SPECTROPORT mobile metal analyzer

### When you can't afford metal mistakes: On-the-spot sorting, identification, verification, and analysis of metal alloys

- Ideal for on-site measurement of metal composition
- Highly cost-effective versus prohibitive expenses of off-spec products
- Pass/fail sorting mode, to separate samples of a given grade or a specific heat
- Grade identification mode, to determine a grade and ensure sample meets specified limits
- Grade verification/control mode, to confirm sample fits preselected grade specifications
- Unique iCAL one-sample standardization, to save average of 30 minutes per day

When it comes to metals specifications, the news is full of mismeasurements, mistakes, and scandals in a host of industries. It's increasingly clear that quality-conscious organizations can't afford to give off their responsibility for metals verification. Inspecting the metals makeup of incoming and outgoing components has become a critical quality control (QC) task for companies worldwide.



Fortunately, inspection can be easily, accurately, and affordably accomplished. Use advanced solutions from SPECTRO Analytical Instruments: either the flagship **SPECTROTEST mobile metal analyzer**, or the more portable **SPECTROPORT metal analyzer**.

Alloys must be accurately sorted, identified, and verified to meet the chemical composition specified by a customer or an industry. For organizations performing metal producing, processing, recycling, or service contracting, an alloy mixup at the shipping dock or on the factory floor risks an expensive, inconvenient batch rework — or a catastrophic lost contract.

Continuous QC can help establish the identity of each metal or alloy, from initial melt to finished product or final application. The process is simple:

To perform *rapid pass/fail sorting or grade identification* with either SPECTRO analyzer, choose arc mode, hold the probe to a sample, and push the start button. Results appear within 2 seconds.

For more demanding *grade verification or analysis* — including detection of additional elements such as carbon, phosphorus, sulfur, and more — select spark mode, so the probe head is purged with argon gas. Again, hold the probe to a sample and push start. Results typically displayed within 10 seconds.



#### GERMANY


SPECTRO Analytical Instruments GmbH  
Boschstrasse 10  
D-47533 Kleve  
Tel: +49.2821.892-0  
Fax: +49.2821.8922200  
spectro.sales@ametec.com

#### U.S.A.

SPECTRO Analytical Instruments Inc.  
91 McKee Drive  
Mahwah, NJ 07430  
Tel: +1.800.548.5809  
+1.201.642.3000  
Fax: +1.201.642.3091  
spectro-usa.sales@ametec.com

## SPECTROTEST & SPECTROPORT mobile metal analyzer

### When you can't afford metal mistakes: On-the-spot sorting, identification, verification, and analysis of metal alloys



SPECTROPORT: Simple and rapid Pass/Fail sorting on a micro alloyed steel.


Spark Analyzer Pro Mobile - Routine Measurement - V 1.06.007 - [Cont.]


Sample	Grade Norm	Grade ID
Stainless steel	DIN	1.4571


	C	Si	Mn	P	S	Cr	Mo	Ni	Al	Co	Cu
1	0.0199	0.554	1.96	0.0335	0.0239	16.67	2.10	11.04	0.0114	0.178	0.400
2	0.0190	0.557	1.96	0.0319	0.0247	16.70	2.09	11.05	0.0122	0.178	0.394
3	0.0199	0.556	1.96	0.0339	0.0239	16.73	2.09	11.04	0.0123	0.178	0.395
IC						16.50	2.00	10.50			
IC	0.0196	0.555	1.97	0.0331	0.0242	16.70	2.09	11.05	0.0120	0.178	0.396
ID	0.0800	1.000	2.00	0.0450	0.0300	18.50	2.50	13.50			


	Nb	Ti	V	W	Pb	Sn	As	Zr	Ca	B	Fe
1	0.0161	0.146	0.0596	<0.0400	<0.0100	0.0159	0.0067	0.0032	0.0011	0.0028	66.7
2	0.0155	0.144	0.0593	<0.0400	<0.0100	0.0144	0.0078	0.0032	0.0009	0.0028	66.7
3	0.0147	0.147	0.0592	<0.0400	<0.0100	0.0153	0.0070	0.0032	0.0010	0.0027	66.6
IC			0.0979								
IC	0.0155	0.146	0.0594	<0.0400	<0.0100	0.0152	0.0071	0.0032	0.0010	0.0028	66.7
ID		0.800									


Grade Norm	Grade Name	Grade Description
DIN	1.4571	X 6 CrNiMoTi 17 12 2


 Stat Values (R)


 iCAL (F)


 Flush (F)


 Recall (S)


 Print (F6)


 Meas Step (F3)


 Export (F5)

 Apply Grade (A)






 Sich Grade (F4)

 Mean/RSD (V)

 Help

 Back

Fe-30 Mo Fe-0Mo | Quantitative Analysis

Ready     

SPECTROTEST: Grade verification in spark mode on a stainless steel sample.

Both instruments routinely perform analyses for ferrous, aluminum, copper, nickel, cobalt and titanium alloys. SPECTROTEST's no-compromise, high-resolution optical system can handle all elements necessary for a complete metal analysis on the spot. So it also delivers excellent results when measuring nitrogen in ferrous-based alloys such as duplex steels; detecting small amounts of lithium and sodium in aluminum-based materials; or it is used also for tin, zinc, lead and magnesium based alloys.

Both analyzers allow operators to change modes quickly and easily, without tools. Adapters are available for the analysis of tubes, wires, and small parts, as well as for special pass/fail forms. SPECTROPORT and SPECTROTEST can operate independently from AC line power. An optional rechargeable battery pack handles up to 800 measurements.

You can't afford a metal mistake. Make sure it's right, every time, with a SPECTRO metal analyzer.